



## Masy Folding Validation Cart



Make moving the full Kaye Validator 2000 package — Validator and laptop computer — easy and effortless when you rent or purchase a folding validation cart from Masy BioServices.

At only 4.9 pounds and 8" wide when collapsed, Masy's folding validation cart is the perfect add-on accessory for your Validator package rental. Easy to unpack and open, the cart locks open for secure, stable moving of your Validator set-up from project to project. Three sturdy shelves hold the laptop, Validator/SIMs thermocouple harness and uninterruptible power supply at the proper working height for efficient chamber mapping.

Constructed of durable, lightweight aluminum with 3-inch rubber swivel wheels for smooth movement.

Maximum load:	65 lbs. per level; 195 lbs. total.
Open Dimensions:	34" X 15" X 26-1/8"
Closed Dimensions:	42" X 15" X 8"

## EQUIPMENT RENTAL & SALES

- ▶ Validation Datalogging Systems
  - GE Kaye Validator 2000
  - DataTrace Wireless Dataloggers
- ▶ Temperature Standards
  - GE Kaye IRTD-400
  - Fluke 1523/1524
- ▶ Temperature Reference Baths
  - GE Kaye LTR-25/140
  - GE Kaye LTR-40/140
  - GE Kaye LTR-90
  - Fluke 9190
  - GE Kaye HTR-400
  - GE Kaye CTR-80
  - Hart 7380
- ▶ Specialty Reference Baths
  - GE Kaye K140-4
  - Hart Micro-Bath 6102
  - Hart Micro-Bath 7103
  - Hart Portable Furnace 9150
  - Hart Handheld Drywell 9100
  - Hart Handheld Drywell 9102

### *Additional Equipment*

- ▶ Rosemount Pressure Transducer
- ▶ International Light UV Wide-Band A
- ▶ International Light Illuminance
- ▶ Monarch Tach 4 Tachometer
- ▶ Multimeter
- ▶ GEO 2000 Humidity Generator
- ▶ Vaisala HMI31 RH/Temperature Meter
- ▶ Hygroclip RH Assembly / HCS + HC2
- ▶ Hygrologger Temperature/RH Assembly
- ▶ Vaportron HL-100CL
- ▶ Ultra-Premium Thermocouples
- ▶ Masy Rev K Thermocouple Feed-Through
- ▶ Masy Rev M Thermocouple Feed-Through
- ▶ Steam Quality Kit
- ▶ RubberFab Gaskets & Clamps
- ▶ CO<sub>2</sub> Sensor Assembly
- ▶ Folding Validation Cart

# Masy Folding Validation Cart

To unfold the cart, pull the handles downwards and outwards, then push down on the brace to lock into place. The cart collapses easily by lifting up on the locking brace and folding the sides upwards and inwards.



## Complete Validator 2000 Rental Packages Include:

### Kaye Validator 2000 Thermal Mapping System with Laptop

The Kaye Validator 2000® is a stand-alone thermal validation system for life science industries – pharmaceuticals, medical devices, and biotechnology – and other regulated and non-regulated applications where precise thermal mapping can affect the quality of product development and storage.



### Kaye IRTD 400

The intelligent resistant temperature device (IRTD) is a secondary standard with a completely self-contained measurement system, including electronics for calibration and temperature conversion. Kaye's IRTD communicates accurate temperature data directly to the Kaye Validator with no additional measurement or readout equipment required.



### Fluke Metrology Wells

The Fluke Calibration 9190A Ultra-Cool Field Metrology Well is the most accurate and stable, cold temperature dry-block on the market.



### Validation Thermocouple Assemblies

Masy BioServices' thermocouples are specifically designed for validation and thermal mapping. We offer thermocouples made from ultra-premium select grade wire; custom pre-assembled harnesses and raw wire spools are available.



Contact Masy BioServices today for a full Kaye Validator 2000 thermal mapping package, including custom-harnessed ultra-premium thermocouples and a convenient folding cart.